
Installing a Robot onto a Rail

Overview


- 1. Choose your configuration
- 2. Route cables
- 3. Attach robot and cables
- 4. Install new PAC files to include linear rail
- 5. Calibrate robot
- 6. Set chosen configuration in PAC files

1. Choose your configuration

O degrees: Robot's Y axis aligns with rail travel




-90 degrees: Robot's X axis aligns with rail travel



2. Route Cables

 If necessary, remove the covers on the top of the linear axis carriage and re-route the cabling so that the connectors can reach the facilities panel at the base of the z-axis column.



3. Attach robot and cables

- 1. Align outer holes of robot with rail
- 2. Attach cables
 - 1. Power
 - 2. Ethernet
 - 3. Communication (DB9)




0 deg: base aligned with edge



-90 deg: base centered

4. Install new PAC files

- 1. Follow standard procedure
 - 1. Make sure to match the version number of the new files to the files currently installed on the robot!



- 2. Reboot robot after new PAC files are installed
- 3. Check Virtual Pendant to see if the rail now shows position

5. Calibrate the robot

- 1. Follow standard procedure
- 2. Rail calibration position is described in the Cal_pp dialogue

6. Set configuration in PAC files

- 1. Set the 5th value of parameter 16050 to the chosen configuration (0 or -90)
- Click "Set new values"
- Click "Save All to Flash"
- 4. Wait 20 seconds
- 5. Reboot the robot

The rail is now ready to use

